


## Camera Tube Cover



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	CTC-001	Camera tube cover	200pcs

### Features:

- Overall size (250x13cm)
- telescopically folded (13x50cm)
- Perforated, tear-off tip
- Inclusive of adhesive tape (2.5x25cm)