


Cardiovascular Pack 1



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	CVP-001a	Contents:	2packs
		1 Reinforced gowns XL with OT Towel in outer wrap	
		2 OT towel	
		2 Reinforced gowns XL	
		1 Reinforced overhead table cover (200x220cm)	
		1 Suture bag	
		6 Small adhesive drape (40x65cm)	
		1 Peri towel adhesive (19x55cm)	
		1 Plastic sheet (120x150cm)	
		1 Bottom sheet absorbent reinforced (137x165cm)	
		2 Side sheet absorbent reinforced, adhesive (100x196cm)	
		1 Head sheet absorbent reinforced, adhesive (132x251cm)	
		1 Foot sheet absorbent reinforced,	

adhesive (132x191cm)
1 Table cover (120x200cm)

Drape features:

- Absorbent impervious fabric
 - Head (64x36cm) - Side (30x185cm) - Foot (63x36cm)
- Cord/tube holders on head and side sheets