


Cardiovascular Pack 2



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	CVP-002	Contents:	4packs
		1 Mayo stand cover (75x145cm)	
		4 OT towel (50x50)	
		2 Foot cover (25x43cm)	
		1 Diathermy pouch	
		2 OR fabric tape (10x50cm)	
		1 Plain drape (150x180cm)	
		1 Split sheet with perineal drape (200x260cm)	
		1 Cardiovascular drape (300x300cm)	
		1 Table cover (150x180cm)	

Drape features:

- Absorbent impervious fabric (120x198cm)

- Full incise (31x40cm)
- Diathermy pouch (34x38cm)
- Armboard covers(31x40cm)
- Cable holder(7x40cm)
- 2 side pouches(compartmented)
- Self-adhesive U-slit perineal drape