


Cesarean Pack with 360 Degree Pouch



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	LCP-003R	Contents:	8 packs
		1 Mayo stand cover (58x140cm)	
		1 Suture Bag	
		2 OT towels (50x50cm)	
		1 LSCS Drape with 360 pouch	
		1 Table cover (120x200cm)	

Drape features:

- Absorbent impervious fabric (43x56cm)
- Incise area (30x30cm) with pear shaped fenestration (20x20cm)
- 360 fluid collection pouch with malleable wires
- Cord/tube holders

