


Laparoscopy / Pelviscopy Pack



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	LPP-001	Contents:	8 packs
		1 Mayo stand cover (58x140cm)	
		1 Suture Bag	
		1 Small adhesive drape (40x65cm)	
		1 Impervious under buttock drape (90x125cm)	
		2 Leggings (72x120cm)	
		1 Laparoscopy / Pelviscopy drape	
		1 Table cover (120x200cm)	

Drape features:

- Absorbent impervious fabric (35x70cm)
- Adhesive trapezium fenestration (17x18x23.5cm)
- Perineal fenestration (8x15cm)

- Attached leggings
- Cord/tube holders