


## Nephroscopy Pack



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	PCNLP-001	Contents: 1 Suture bag 1 OT Towel (50x50cm) 1 Nephroscopy drape with pouch 1 Table cover (120x200cm)	8 packs

### Drape features:

- Absorbent impervious fabric (90x120cm)
- Adhesive fenestration (19x25cm)
- Fluid collection pouch (61x100cm)
- Cord/tube holders