


Arthroscopy Drape



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	ATD-004	Arthroscopy Drape	12pcs

Drape features:

- Absorbent impervious fabric (50x57cm)
- 6cm (Top), 5cm (bottom) elastic fenestration
- Attached with tying string fluid collection pouch
- Cord/tube holders