


Cardio-Thoracic Drape



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description Cat. Code	Description	Qty. / Ctn.
CTD-001	Cardio-Thoracic Drape	8pcs

**Drape features:**

- Absorbent impervious fabric (120x198cm)
- Full incise (31x40cm)
- Diathermy pouch (34x38cm)
- Armboard covers
- Cable holder (7x40cm)
- 2 sided pouches(compartmented)