


Cesarean Drape with 270 Degree Pouch 2



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description Cat. Code	Description	Qty. / Ctn.
LSD-003U	Cesarean drape with 270° pouch incise fenestration	7pcs

Drape features:

- Absorbent impervious fabric (43x56cm)
- Incise area (30x30cm) with pear shaped fenestration (20x20cm)
- 270° fluid collection pouch with malleable wires
- Cord/tube holders