


Laparoscopic Assisted Vaginal Hysterectomy Drape



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	LVD-001a	Laparoscopic assisted vaginal hystrectomy	10pcs

Drape features:

- Absorbent impervious fabric (35x92cm)
- Adhesive trapezium fenestration (17x23x23cm)
- Perineal adhesive fenestration (8x15cm)
- Instrument pad
- Attached with tying string fluid collection pouch
- Attached leggings
- Armboard covers
- Flexible cord/tube holders