


Nephroscopy Drape with Pouch



Rating: Not Rated Yet

Price:

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	PCNLP-001	Nephroscopy drape with pouch	12pcs

Drape features:

- Absorbent impervious fabric (90x120cm)
- Adhesive fenestration (19x25cm)
- Fluid collection pouch (61x100cm)
- Cord/tube holders