


## Transuretral Drape with Pouch



Rating: Not Rated Yet

**Price:**

Variant price modifier:

Base price with tax:

Price with discount:

Salesprice with discount:

Sales price:

Sales price without tax:

Discount:

Tax amount:

[Ask a question about this product](#)

Manufacturer: [hopes](#)

Description	Cat. Code	Description	Qty. / Ctn.
	TUR-001	Transurethral drape with pouch	10pcs

### Drape features:

- Absorbent impervious fabric (60x78cm)
- Adhesive abdominal fenestration (8cm ?)
- Perineal aperture (3cm ?)
- Finger cot
- Fluid-control pouch with filtering mesh